

Technical Datasheet

		EN ISO 13485 certified
Article - No:	04-096-6008	Revision: 7
Remarks		single use only

Specification

Product	Adhesive sealing film, for cold-storage, 6-strip, 50µm, -80°C up to +120°C, alu., pierceable, np pcr ready, sterile R
Dimensions	length: 101,00 mm (+/- 0,2 mm)
	width: 60,00 mm (+/- 0,2 mm) – 1 strip = 10,00 mm
Material Characteristics	aluminium (AL)
Temperature range	-80°C - 120°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 300
	pcs per carton: 300
	pcs per pallet: 700.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department