

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-812-8062	Revision: 1
Remarks		single use only

Specification

Product	PS stick with big viscose tip, Ø2,5x150mm, single peel-pack, np pcr ready, sterile EO, CE2797
Classification	CE2797
Dimensions	diameter: 2,50 mm (+/- 0,2 mm)
	length: 150,00 mm (+/- 0,2 mm)
	diameter tip : 4,00-5,50 mm (+/- 0,2 mm)
Material Characteristics	stick: polystyrene (PS)
	tip: viscose (VI)
Product Characteristics	biologically inert
	physiologically harmless
	free of detectable phthalate acc. to REACH regulation
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	gas sterilisation by ethyleneoxid
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 100
	pcs per carton: 1.200
	pcs per pallet: 128.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

