

Technical Datasheet

Date: 22.06.20

Page 1 of 1

		EN ISO 13485 certified
Article - No:	09-819-0102	Revision: 2
Remarks		single use only

Specification

Product	Wooden applicator with small cotton tip, Ø2,2x150mm, np pcr ready, CE
Classification	CE
Dimensions	diameter stick: 2,20 mm (+/- 0,2 mm)
	length stick: 150,00 mm (+/- 0,2 mm)
Material Characteristics	stick: wood (TI)
	tip: cotton (CO)
Product Characteristics	biologically inert
	physiologically harmless
	free of detectable phthalate acc. to REACH regulation
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	N/A
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 100
	pcs per carton: 10.000
	pcs per pallet: 560.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	