

SENSIColumn Immunoaffinity Columns



Product	Product description	Quantity	Article No.
SENSIColumn IAC Aflatoxin total	Immunoaffinity columns that retain total aflatoxin	3 ml x 50 columns	EFAF3105
SENSIColumn IAC Aflatoxin M1	Immunoaffinity columns that retain aflatoxin M1	3 ml x 50 columns	EFAFM3115
SENSIColumn IAC Ochratoxin	Clean-up of commodity extracts of food and feed samples containing ochratoxin A/B	3 ml x 50 columns	EFOT3125
SENSIColumn IAC Ochratoxin SMART	Specially designed low cost immunoaffinity column with reduces loding capacity for ochratoxin purification.	3 ml x 50 columns	EFOTS3125
SENSIColumn IAC Deoxynivalenol	Immunoaffinity columns that retain deoxynivalenol and separate it from the other substances present in the sample.	3 ml x 50 columns	EFD03135
SENSIColumn IAC Fumonisin	Immunoaffinity columns that retain fumonisins and separate them from the other substances present in the sample.	3 ml x 50 columns	EFFU3145
SENSIColumn IAC Zearalenone	Immunoaffinity columns that retain zearalenone and separate it from the other substances present in the sample.	3 ml x 50 columns	EFZE3155
SENSIColumn IAC Zearalenone SMART	Specially designed low cost immunoaffinity column with reduces loding capacity for ochratoxin purification.	3 ml x 50 columns	EFZES3155
SENSIColumn IAC T2/HT2	Immunoaffinity columns that retain T2 toxin and separate it from the other substances present in the sample.	3 ml x 50 columns	EFTE3165
SENSIColumn IAC Combi Afla+Ochra	Clean-up of commodity extracts of food and feed samples containing aflatoxins and ochratoxins	3 ml x 50 columns	EFCA3215
SENSIColumn IAC Combi Afla+Ochra+Zeara	Clean-up commodity extracts of food and feed samples containing aflatoxins, ochratoxins and zearalenon	3 ml x 50 columns	EFCZ3225
SENSIColumn IAC Combi Deoxinivalenol+Zearalenone	Clean-up of commodity extracts of food and feed samples containing deoxynivalenol and zearalenon	3 ml x 50 columns	EFGD3235
SENSIColumn IAC Combi AOZDFT all in one for LC-MS/MS	Simultaneous extraction of aflatoxins, ochratoxins, zearalenon, deoxynivalenol, fuminisins and T2-Toxin in food and feed stuffs	3 ml x 50 columns	EFCM3265
SENSIColumn IAC Ergot Alkaloids	Clean-up of commodity extracts of food and feed samples containing 12 priority ergot alkaloids	3 ml x 50 columns	EFLS3275

Distributed by



info@neotec srl 0817599062